


Product Datasheet

MAT2B antibody (orb12330)

Description	Goat polyclonal to MAT2B.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MAT2B (isoform 1)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize reported isoform 1 NP_037415.1 only.
Protein Sequence	REKELSIHFVPGS
Clonality	Polyclonal
MW	37.6; 36.4
Entrez	27430
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Thymus ly...