

Product Datasheet

USP6 Peptide (orb1233393)

Description	USP6 Peptide
Tested Applications	Blocking
Target	USP6
Preservatives	USP6 peptide is supplied as 100 µg of dried peptide. When the peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL.
Form/Appearance	Lyophilized
Storage	Store at -20°C. When reconstituted minimize freezing and thawing.
Note	For research use only
Application notes	USP6 Peptide is used for blocking the activity of USP6 antibody (catalog number 46-565). The USP6 antibody ships at 0.5 mg/mL and the USP6 Peptide ships as a 100 µg pellet. When peptide is reconstituted in 200 µL water, the concentration will be 0.5 mg/mL, identical to that of the antibody. The best starting ratio would be 1:1 (which means molar excess of peptide relative to antibody when identical volumes are mixed). Mix equal volume of peptide and antibody at the required dilution and leave at ambient temperature for 20 minutes. Then incubate and develop in parallel two blots: one with the antibody and one with the antibody which has been pre-absorbed with the provided peptide.
Uniprot ID	P35125
NCBI	NP_004496
Expiration Date	6 months from date of receipt.