

Product Datasheet

Ryk antibody (orb12336)

Description	Goat polyclonal to Ryk.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Ryk (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_038677.3;NP_001036072.1.
Protein Sequence	HAKSKVEYKLGFG
Clonality	Polyclonal
MW	66.3; 66.1
Entrez	6259
Dilution Range	ELISA: 1:32000, WB: 0.1 µg/ml
Expiration Date	12 months from date of receipt.