

Product Datasheet

ATG10 Peptide (orb1234721)

Description	ATG10 Peptide
Tested Applications	Blocking
Target	ATG10
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store ATG10 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	ATG10 peptide is used for blocking the activity of ATG10 antibody.
Uniprot ID	Q9H0Y0
NCBI	EAW95884
Expiration Date	6 months from date of receipt.