


Product Datasheet

Cyb5r4 antibody (orb12356)

| | |
|----------------------------|--|
| Description | Goat polyclonal to Cyb5r4. |
| Species/Host | Goat |
| Reactivity | Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Target | Cyb5r4 (mouse aa444-458) |
| Preservatives | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing. |
| Storage | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Note | For research use only |
| Application notes | WB: Approx 65kDa band observed in Rat Spleen lysates (calculated MW of 55kDa according to NP_596918.2 and 58.8kDa according to Q68EJ0-1). Recommended concentration: 0.1-0.3µg/ml. |
| Protein Sequence | RCQLEKLALREKRFD |
| Clonality | Polyclonal |
| MW | 59.7 |
| Dilution Range | ELISA: 1:64000, WB: 0.1-0.3 µg/ml |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Rat Spleen lysa...