

Product Datasheet

UEV1A Peptide (orb1236658)

Description	UEV1A Peptide
Tested Applications	Blocking
Target	UBE2V1
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store UEV1A peptide at -20°C, stable for one year.
Note	For research use only
Application notes	UEV1A peptide is used for blocking the activity of UEV1A antibody.
Uniprot ID	Q13404
NCBI	NP_954595
Expiration Date	6 months from date of receipt.