

Product Datasheet

CARMA3 Peptide (orb1236712)

Description	CARMA3 Peptide
Tested Applications	Blocking
Target	CARD10
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store CARMA3 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	CARMA3 peptide is used for blocking the activity of CARMA3 antibody.
Uniprot ID	Q9BWT7
NCBI	NP_055365
Expiration Date	6 months from date of receipt.