

Product Datasheet

Caspase-13 Peptide (orb1236825)

Description	Caspase-13 Peptide
Tested Applications	Blocking
Target	CASP4
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store Caspase-13 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	Caspase-13 peptide is used for blocking the activity of caspase-13 antibody.
Uniprot ID	O75601
NCBI	O75601
Expiration Date	6 months from date of receipt.