

## Product Datasheet

# Goat Anti-Rabbit IgG Antibody (H&L) Antibody (HRP) (orb1236959)

<b>Description</b>	Goat Anti-Rabbit IgG Antibody (H&L) Antibody (HRP)
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, ICC, IHC, WB
<b>Immunogen</b>	Purified rabbit IgG.
<b>Preservatives</b>	Buffer containing 0.01M sodium phosphate, 0.25M sodium chloride, pH 7.6 and 25% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store 2 - 8 °C as an undiluted liquid for several weeks. After dilution, do not use for more than one day. -20 °C for up to one year.
<b>Note</b>	For research use only
<b>Application notes</b>	1:250-1:3,000 for IHC/ICC1:2,000-1:5,000 for ELISA, Western Blot, color substrates.1:5,000-1:10,000 for Western Blot, chemiluminescence.
<b>Isotype</b>	IgG
<b>Dilution Range</b>	1:250-1:3,000 for IHC/ICC1:2,000-1:5,000 for ELISA, Western Blot, color substrates.1:5,000-1:10,000 for Western Blot, chemiluminescence.
<b>Expiration Date</b>	12 months from date of receipt.