
Product Datasheet

**Mouse Anti-Rabbit IgG Antibody [6C1A8]
(orb1236971)**

Description

Mouse Anti-Rabbit IgG Antibody [6C1A8]

Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Rabbit IgG
Preservatives	PBS pH 7.2, 1% rAlbumin, 0.02% Thimerosal
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 2-8°C. The product is stable for about 4-6 weeks at 4°C as an undiluted liquid. For extended storage, aliquot and freeze the product at -70°C or below. Avoid repeated freezing and thawing. Or alternatively, add an equal volume of glycerol for a final concentration of 50% and store at -20°C.
Note	For research use only
Application notes	Immunoblotting: 1/5,000 - 1/50,000, ELISA: 1/5,000 - 1/50,000, Immunohistochemistry: 1/100 - 1/1,000
Isotype	IgG2
Dilution Range	Immunoblotting: 1/5,000 - 1/50,000, ELISA: 1/5,000 - 1/50,000, Immunohistochemistry: 1/100 - 1/1,000
Expiration Date	12 months from date of receipt.