
Product Datasheet

**Mouse Anti-Goat IgG Antibody [1G12H5] (biotin)
(orb1236980)**

Description

Mouse Anti-Goat IgG Antibody

Species/Host	Mouse
Reactivity	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Goat IgG
Preservatives	PBS pH 7.2, 1% rAlbumin, 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 2-8°C. The product is stable for about 4-6 weeks at 4°C as an undiluted liquid. For extended storage, aliquot and freeze the product at -70°C or below. Avoid repeated freezing and thawing. Or alternatively, add an equal volume of glycerol for a final concentration of 50% and store at -20°C.
Note	For research use only
Application notes	Immunoblotting: 1/1,000 - 1/20,000, ELISA: 1/1,000 - 1/20,000, Immunohistochemistry: 1/100 - 1/1,000
Isotype	IgG2
Dilution Range	Immunoblotting: 1/1,000 - 1/20,000, ELISA: 1/1,000 - 1/20,000, Immunohistochemistry: 1/100 - 1/1,000
Expiration Date	12 months from date of receipt.