

---

## Product Datasheet

**Mouse Anti-Rabbit IgG Antibody [6C1A8] (HRP)  
(orb1236983)**

**Description**

Mouse Anti-Rabbit IgG Antibody

<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Rabbit IgG
<b>Preservatives</b>	PBS pH 7.2, 1% rAlbumin, 0.02% Thimerosal
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Store at 2-8°C. The product is stable for about 4-6 weeks at 4°C as an undiluted liquid. For extended storage, aliquot and freeze the product at -70°C or below. Avoid repeated freezing and thawing. Or alternatively, add an equal volume of glycerol for a final concentration of 50% and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Immunoblotting: 1/5,000 - 1/50,000, ELISA: 1/5,000 - 1/50,000, Immunohistochemistry: 1/100 - 1/1,000
<b>Isotype</b>	IgG2
<b>Dilution Range</b>	Immunoblotting: 1/5,000 - 1/50,000, ELISA: 1/5,000 - 1/50,000, Immunohistochemistry: 1/100 - 1/1,000
<b>Expiration Date</b>	12 months from date of receipt.