

## Product Datasheet

# Control Normal Hamster IgG (APC-Cy5.5) Antibody (orb1237039)

<b>Description</b>	Control Normal Hamster IgG (APC-Cy5.5) Antibody
<b>Species/Host</b>	Hamster
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Preservatives</b>	Supplied in PBS/sodium azide and a stabilizing agent
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.1 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Normal Hamster IgG (APC-Cy5.5) can be used as control in flow cytometry assays.
<b>Isotype</b>	IgG
<b>Dilution Range</b>	Normal Hamster IgG (APC-Cy5.5) can be used as control in flow cytometry assays.
<b>Expiration Date</b>	12 months from date of receipt.