

## Product Datasheet

# Control Normal Hamster IgG (APC) Antibody (orb1237044)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Control Normal Hamster IgG (APC) Antibody                                 |
| <b>Species/Host</b>        | Hamster   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | FC, FLISA   |
| <b>Preservatives</b>       | Supplied in PBS/sodium azide and a stabilizing agent                      |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 0.1 mg/mL   |
| <b>Storage</b>             | Store vial at 2-8°C   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Normal Hamster IgG (APC) can be used as control in flow cytometry assays. |
| <b>Isotype</b>             | IgG   |
| <b>Dilution Range</b>      | Normal Hamster IgG (APC) can be used as control in flow cytometry assays. |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |