

Product Datasheet

Control Normal Goat F(ab')₂ IgG (APC) Antibody (orb1237121)

Description	Control Normal Goat F(ab') ₂ IgG (APC) Antibody
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	FC, FLISA
Preservatives	Supplied in PBS/sodium azide and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Normal Goat F(ab') ₂ IgG (APC) can be used as control in flow cytometry assays.
Isotype	IgG
Dilution Range	Normal Goat F(ab') ₂ IgG (APC) can be used as control in flow cytometry assays.
Expiration Date	12 months from date of receipt.