

Product Datasheet

Control Normal Goat IgG (APC) Antibody (orb1237128)

Description	Control Normal Goat IgG (APC) Antibody
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	FC, FLISA
Preservatives	Supplied in PBS/sodium azide and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Normal Goat IgG (APC) can be used as control in flow cytometry assays.
Isotype	IgG
Dilution Range	Normal Goat IgG (APC) can be used as control in flow cytometry assays.
Expiration Date	12 months from date of receipt.