


Product Datasheet

Zanidatamab (ERBB2/EGFR2/CD340) - Research Grade Biosimilar (orb1237245)

Description	Zanidatamab (ERBB2/EGFR2/CD340) - Research Grade Biosimilar
Conjugation	Unconjugated
Immunogen	Humanized / ERBB2/EGFR2/CD340 [Homo sapiens]
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-kappa/scFv-h-CH2-CH3
Clonality	Recombinant
Source	CHO cells
Expiration Date	12 months from date of receipt.



SDS-PAGE Image

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)