


## Product Datasheet

# Urelumab (TNFRSF9/CD137 ) - Research Grade Biosimilar Antibody (orb1237284)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Urelumab (TNFRSF9/CD137 ) - Research Grade Biosimilar Antibody  |
| <b>Conjugation</b>     | Unconjugated  |
| <b>Immunogen</b>       | Homo sapiens / TNFRSF9/CD137 [Homo sapiens]   |
| <b>Preservatives</b>   | PBS buffer pH 7.5   |
| <b>Concentration</b>   | batch dependent   |
| <b>Storage</b>         | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage. |
| <b>Note</b>            | For research use only   |
| <b>Isotype</b>         | IgG4-kappa  |
| <b>Clonality</b>       | Recombinant   |
| <b>Source</b>          | CHO cells   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |



Urelumab binds with 4-1BB



SDS-PAGE Image