


Product Datasheet

Tucotuzumab Celmoleukin (EPCAM/CD326) - Research Grade Biosimilar Antibody (orb1237288)

| | |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Tucotuzumab Celmoleukin (EPCAM/CD326) - Research Grade Biosimilar Antibody |
| Conjugation | Unconjugated |
| Immunogen | Humanized / EPCAM/CD326 [Homo sapiens] |
| Preservatives | PBS buffer pH 7.5 |
| Concentration | batch dependent |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage. |
| Note | For research use only |
| Isotype | IgG1-nd |
| Clonality | Recombinant |
| Source | CHO cells |
| Expiration Date | 12 months from date of receipt. |



SDS-
PAGE
Image