

Product Datasheet

Solitomab (CD3E/EpCAM) - Research Grade Biosimilar Antibody (orb1237358)

Description	Solitomab (CD3E/EpCAM) - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Mus musculus / CD3E/EpCAM [Homo sapiens]
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Isotype	(scFv-kappa-heavy)-(scFv-heavy-kappa)
Clonality	Recombinant
Source	CHO cells
Expiration Date	12 months from date of receipt.