


Product Datasheet

Remtolumab (IL17A/TNFSF2/TNF-alpha/TNFA) - Research Grade Biosimilar Antibody (orb1237410)

Description	Remtolumab (IL17A/TNFSF2/TNF-alpha/TNFA) - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Homo sapiens / IL17A/TNFSF2/TNF-alpha/TNFA [Homo sapiens]
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Isotype	[VH-VH'-H-Gamma1-VL-VL'-C-kappa]-dimer
Clonality	Recombinant
Source	CHO cells
Expiration Date	12 months from date of receipt.



SDS-PAGE Image