


Product Datasheet

Racotumomab (Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) - gangliosides GM3) Mus musculus IgM-kappa monoclonal antibody P3) - Research Grade Biosimilar (orb1237424)

Description	Racotumomab (Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) - gangliosides GM3) Mus musculu...
Conjugation	Unconjugated
Immunogen	Mus musculus / Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) - gangliosides GM3) Mus musculus IgM-kappa monoclonal antibody P3 [Homo sapiens]
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Isotype	IgG1-kappa
Clonality	Recombinant
Source	CHO cells
Expiration Date	12 months from date of receipt.



SDS-PAGE Image