


Product Datasheet

Iodine (131I) Derlotuximab Biotin (DNA/histone(H1) complex) - Research Grade Biosimilar Antibody (orb1237610)

| | |
|------------------------|---|
| Description | Iodine (131I) Derlotuximab Biotin (DNA/histone(H1) complex) - Research Grade Biosimilar Antibody |
| Conjugation | Unconjugated |
| Immunogen | Chimeric / DNA/histone(H1) complex [Homo sapiens] |
| Preservatives | PBS buffer pH 7.5 |
| Concentration | batch dependent |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage. |
| Note | For research use only |
| Isotype | IgG1-kappa |
| Clonality | Recombinant |
| Source | CHO cells |
| Expiration Date | 12 months from date of receipt. |



SDS-PAGE Image