

## Product Datasheet

# Daratumumab (CD38) - Research Grade Biosimilar Antibody (orb1237758)

**Description** Daratumumab (CD38) - Research Grade Biosimilar Antibody

**Conjugation** Unconjugated

**Immunogen** Homo sapiens / CD38 [Homo sapiens]

**Preservatives** PBS buffer pH 7.5

**Concentration** batch dependent

**Storage** Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.


**Note** For research use only

**Isotype** IgG1-kappa


**Clonality** Recombinant

**Source** CHO cells

**Expiration Date** 12 months from date of receipt.



Daratumumab binds with CD38



SDS-PAGE Image