

Product Datasheet

CD3E Biosimilar Antibody (orb1237909)

Description	CD3E Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Chimeric / CD3E, CD33
Target	CD3E
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific Homodimer (VK-VH-VL'-VH', Tandem diabody);Kappa;Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2243775-32-6
Uniprot ID	P07766
NCBI	P07766
Expiration Date	12 months from date of receipt.