

Product Datasheet

C5 Biosimilar Antibody (orb1237911)

| | |
|------------------------|--|
| Description | C5 Biosimilar Antibody |
| Conjugation | Unconjugated |
| Immunogen | Chimeric / C5 |
| Target | C5 |
| Preservatives | PBS buffer pH 7.5 |
| Concentration | batch dependent |
| Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term. |
| Note | For research use only |
| Isotype | IgG4,Kappa |
| Clonality | Recombinant |
| Purity | > 95% |
| Source | CHO cells |
| CAS Number | 2250440-41-4 |
| Uniprot ID | P01031 |
| NCBI | P01031 |
| Expiration Date | 12 months from date of receipt. |