

Product Datasheet

Ceacam5 Biosimilar Antibody (orb1237918)

Description	Ceacam5 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Chimeric / CEACAM5
Target	Ceacam5
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2349294-95-5
Uniprot ID	Q3UUK2
NCBI	Q3UUK2
Expiration Date	12 months from date of receipt.