

Product Datasheet

IL31 Biosimilar Antibody (orb1237922)

Description	IL31 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Caninized / IL31 (Canine)
Target	IL31
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG2,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2364504-80-1
Uniprot ID	C7G0W1
NCBI	C7G0W1
Expiration Date	12 months from date of receipt.