

Product Datasheet

PDCD1 Biosimilar Antibody (orb1237926)

Description	PDCD1 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / PDCD1, LAG3
Target	PDCD1
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific scFv ,IgG4;Kappa;Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2245725-04-4
Uniprot ID	Q15116
NCBI	Q15116
Expiration Date	12 months from date of receipt.