

Product Datasheet

pd1 Biosimilar Antibody (orb1237930)

Description	pd1 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Homo sapiens / PDL1
Target	pd1
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG4,Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2256084-03-2
Uniprot ID	B5XDX7
NCBI	B5XDX7
Expiration Date	12 months from date of receipt.