

Product Datasheet

TNFSF13 Biosimilar Antibody (orb1237934)

Description	TNFSF13 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / TNFSF13
Target	TNFSF13
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG2,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2382896-07-1
Uniprot ID	O75888
NCBI	O75888
Expiration Date	12 months from date of receipt.