

Product Datasheet

Rosopatamab - Research Grade Biosimilar Antibody (orb1237937)

Description	Rosopatamab - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Chimeric / FOLH1
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2260767-49-3
Uniprot ID	E9PKM3
NCBI	E9PKM3
Expiration Date	12 months from date of receipt.