

Product Datasheet

CSF2 Biosimilar Antibody (orb1237945)

Description	CSF2 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CSF2
Target	CSF2
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2377482-36-3
Uniprot ID	P04141
NCBI	P04141
Expiration Date	12 months from date of receipt.