

## Product Datasheet

### CD3E Biosimilar Antibody (orb1237951)

<b>Description</b>	CD3E Biosimilar Antibody
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Humanized / TNFRSF17, CD3E
<b>Target</b>	CD3E
<b>Preservatives</b>	PBS buffer pH 7.5
<b>Concentration</b>	batch dependent
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
<b>Note</b>	For research use only
<b>Isotype</b>	Bispecific scFv,Kappa;Lambda
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 95%
<b>Source</b>	CHO cells
<b>CAS Number</b>	2250292-39-6
<b>Uniprot ID</b>	<b>P07766</b>
<b>NCBI</b>	<b>P07766</b>
<b>Expiration Date</b>	12 months from date of receipt.