

Product Datasheet

CD3E Biosimilar Antibody (orb1237953)

Description	CD3E Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / TNFRSF17, CD3E
Target	CD3E
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific scFv,Kappa;Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2251756-52-0
Uniprot ID	P07766
NCBI	P07766
Expiration Date	12 months from date of receipt.