

## Product Datasheet

# Onfekafusp - Research Grade Biosimilar Antibody (orb1237956)

<b>Description</b>	Onfekafusp - Research Grade Biosimilar Antibody
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Homo sapiens / FN extra domain B
<b>Preservatives</b>	PBS buffer pH 7.5
<b>Concentration</b>	batch dependent
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
<b>Note</b>	For research use only
<b>Isotype</b>	Kappa
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 95%
<b>Source</b>	CHO cells
<b>CAS Number</b>	1957239-88-1
<b>Expiration Date</b>	12 months from date of receipt.