

Product Datasheet

Obrindatamab - Research Grade Biosimilar Antibody (orb1237959)

Description	Obrindatamab - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CD276, CD3E
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific scFv with Domain Crossover,Kappa;Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2069959-72-2
Uniprot ID	H0YKQ9
NCBI	H0YKQ9
Expiration Date	12 months from date of receipt.