

Product Datasheet

Lorukafusp - Research Grade Biosimilar Antibody (orb1237970)

Description	Lorukafusp - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Chimeric / Ganglioside GD2
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1κ
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2131168-99-3
Uniprot ID	Q6RYR7
NCBI	Q6RYR7
Expiration Date	12 months from date of receipt.