

Product Datasheet

TNF Biosimilar Antibody (orb1237972)

Description	TNF Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / TNFA
Target	TNF
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	1688648-13-6
Uniprot ID	P01375
NCBI	P01375
Expiration Date	12 months from date of receipt.