

Product Datasheet

CD47 Biosimilar Antibody (orb1237973)

Description	CD47 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Homo sapiens / CD47
Target	CD47
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG4,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2283356-07-8
Uniprot ID	Q08722
NCBI	Q08722
Expiration Date	12 months from date of receipt.