

## Product Datasheet

### ICOS Biosimilar Antibody (orb1237976)

<b>Description</b>	ICOS Biosimilar Antibody
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Humanized / ICOS, PDCD1
<b>Target</b>	ICOS
<b>Preservatives</b>	PBS buffer pH 7.5
<b>Concentration</b>	batch dependent
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
<b>Note</b>	For research use only
<b>Isotype</b>	Bispecific Mixed mAb and scFv IgG1;na,Kappa;Kappa
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 95%
<b>Source</b>	CHO cells
<b>CAS Number</b>	2329669-80-7
<b>Uniprot ID</b>	<a href="#">Q9Y6W8</a>
<b>NCBI</b>	<a href="#">Q9Y6W8</a>
<b>Expiration Date</b>	12 months from date of receipt.