

Product Datasheet

SLC1A5 Biosimilar Antibody (orb1237979)

Description	SLC1A5 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Homo sapiens / SLC1A5
Target	SLC1A5
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1 κ
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2245205-37-0
Uniprot ID	Q15758
NCBI	Q15758
Expiration Date	12 months from date of receipt.