

Product Datasheet

FCER2 Biosimilar Antibody (orb1237981)

Description	FCER2 Biosimilar Antibody
Conjugation	Unconjugated
Target	FCER2
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1 κ
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	357613-86-6
Uniprot ID	P06734
NCBI	P06734
Expiration Date	12 months from date of receipt.