

Product Datasheet

CTLA4 Biosimilar Antibody (orb1237988)

Description	CTLA4 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CD274, CTLA4
Target	CTLA4
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific Single Domains (VH-VH'-CH),IgG1na;na
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2367013-69-0
Uniprot ID	P16410
NCBI	P16410
Expiration Date	12 months from date of receipt.