

Product Datasheet

ADM Biosimilar Antibody (orb1237989)

Description	ADM Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / ADM
Target	ADM
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1κ
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2305638-98-4
Uniprot ID	P35318
NCBI	P35318
Expiration Date	12 months from date of receipt.