

Product Datasheet

Cadonilimab - Research Grade Biosimilar (orb1238007)

Description	Cadonilimab - Research Grade Biosimilar
Conjugation	Unconjugated
Immunogen	Humanized / PDCD1, CTLA4
Target	PDCD1
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific ,IgG1;Kappa;Lambda
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2394841-59-7
Uniprot ID	Q15116
NCBI	Q15116
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)