

Product Datasheet

ERBB2 Biosimilar Antibody (orb1238012)

Description	ERBB2 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / ERBB2 (Domain II), ERBB2 (Domain IV)
Target	ERBB2
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific ,IgG1;Kappa;Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2367012-88-0
Uniprot ID	P04626
NCBI	P04626
Expiration Date	12 months from date of receipt.