

Product Datasheet

TNFSF4 Biosimilar Antibody (orb1238013)

Description	TNFSF4 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / TNFSF4
Target	TNFSF4
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG4,Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2378692-15-8
Uniprot ID	P23510
NCBI	P23510
Expiration Date	12 months from date of receipt.