

Product Datasheet

CD274 Biosimilar Antibody (orb1238016)

Description	CD274 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CD274, TNFRSF9
Target	CD274
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Bispecific IgG1;Lambda;Kappa
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
CAS Number	2253937-12-9
Uniprot ID	Q9NZQ7
NCBI	Q9NZQ7
Expiration Date	12 months from date of receipt.