

## Product Datasheet

### CD3E Biosimilar Antibody (orb1238017)

|                        |  |
|------------------------|--|
| <b>Description</b>     | CD3E Biosimilar Antibody   |
| <b>Conjugation</b>     | Unconjugated   |
| <b>Immunogen</b>       | Chimeric, Humanized / FOLH1, CD3E  |
| <b>Target</b>          | CD3E   |
| <b>Preservatives</b>   | PBS buffer pH 7.5  |
| <b>Concentration</b>   | batch dependent  |
| <b>Storage</b>         | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term. |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | Bispecific scFv;Kappa;Lambda   |
| <b>Clonality</b>       | Recombinant  |
| <b>Purity</b>          | > 95%  |
| <b>Source</b>          | CHO cells  |
| <b>CAS Number</b>      | 2314491-93-3   |
| <b>Uniprot ID</b>      | <b>P07766</b>  |
| <b>NCBI</b>            | <b>P07766</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |